










<https://alodenti.com/>

09927311262

<https://instagram.com/alodenti.shop>

03134341803

Prosthetic Implant Abutment	Item Description	Degree	D	GH	Code	Photo	
	Straight Abutment (SA) Regular	N/A	4.5	1	SA4510		
				2	SA4520		
				3	SA4530		
				4	SA4540		
				5	SA4550		
				6	SA4560		
		N/A	5.5	1	SA5510		
				2	SA5520		
				3	SA5530		
				4	SA5540		
	For more visit our Website : <a href="https://alodenti.com/">https://alodenti.com/</a>						
	Angle Abutment (AA) Regular	15	4.5	1.5	AA154515		
				2.5	AA154525		
				3.5	AA154535		
				4.5	AA154545		
		15	5.5	1.5	AA155515		
				2.5	AA155525		
				3.5	AA155535		
		25	4.5	1.5	AA254515		
				2.5	AA254525		
				3.5	AA254535		
	For more visit our Website : <a href="https://alodenti.com/">https://alodenti.com/</a>						
	UCLA Plastic Abutment Hex & Non Hex	N/A	4.5	1			
	UCLA CCM Abutment Hex & Non Hex						
Normal torque foe gular Abutment is 30N by screwdriver 1.25							
Pre Milled Abutment							








<https://alodenti.com/>

09927311262

<https://instagram.com/alodenti.shop>

03134341803

Prosthetic Implant Abutment

Item Description	Degree	D	GH	Code	Photo
Impression Open Regular	N/A	4.5	2		
	N/A	5.5	2		
	N/A	6.5	2		
Impression Close Regular	N/A	4.5	2		
	N/A	5.5	2		
Normal torque foe gular Abutment is 10N by screwdriver 1.25					
Analog Regular	N/A	4.5			
For more visit our Website : <a href="https://alodenti.com/">https://alodenti.com/</a>					
Healing	N/A	4.5	2.5		
			3.5		
			4.5		
			5.5		
	N/A	5.5	2.5		
			3.5		
			4.5		
			5.5		
N/A	6.5	3.5			
		4.5			
For more visit our Website : <a href="https://alodenti.com/">https://alodenti.com/</a>					
Healing Dentium Abutment	N/A	M	4.5		
			5.5		
			6.5		
	N/A	L	4.5		
			5.5		
			6.5		














<https://alodenti.com/>

09927311262

<https://instagram.com/alodenti.shop>

03134341803

Overdenture

Item Description	Degree	D	GH	Code	Photo
Ball Abutment	2.5	4.5	1		
			2		
			3		
			4		
Castable attachment					
	Normal castable	Green Red		Only for laboratory	Normal castable
	Micro castable				
Oring					
	Normal size	Black		Standard retention	Slightly elastic Maximum suggested time of duration in the mouth : 12 month
	Micro size				
Locator Cap Locator abutment for korean system regular					
GH 1	GH 2	GH 3	GH 4	GH 5	GH 6
					













<https://alodenti.com/>

09927311262

<https://instagram.com/alodenti.shop>

03134341803

Overdenture	Locator				
	Product		Color	Retention	Ref
		Pink cap	Pink	5 IBS	Pink male cap 5 IBS offers the retention you need for your overdenture attachment system
		Light pink cap	Light pink	3 IBS	Light pink male cap 3 IBS offers the retention you need for your overdenture attachment system
		Blue cap	Blue	1.5 IBS	Blue male cap 1.5 IBS offers the retention you need for your overdenture attachment system
		White spacer	White	processing	Locator black out spacer
		Metal Housing	Silver		
	<b>Classic cap size available</b> <b>Normal ( 1.8 mm / 2.2 mm / 2.5mm )</b> <b>Retentive cap color and retention</b>				
	Product		Color	Retention	Ref
		Hard cap Standard cap	Clear	Normal 1500 g Normal 1300 g	Slightly elastic Maximum suggested time of duration in the mouth : 12 month
		Soft cap	Pink	Normal 900 g	elastic Maximum suggested time of duration in the mouth : 12 month
		Extra soft cap	Yellow	Normal 500 g	Very elastic Maximum suggested time of duration in the mouth : 12 month
		Lab cap	Black	Only for laboratory processing	Processing Cap
		Normal Housing	Silver		Normal size 2.2 & 2.5 mm cap



<https://alodenti.com/>

09927311262

<https://instagram.com/alodenti.shop>

03134341803

**Instrument**

**Product**

**Photo**

Analog holder



ball Abutment Screwdriver



Drill Extention



Extracting & insertion tool for removing or inserting the retention caps in the ball Abutment



Extracting & insertion tool for removing or inserting the retention caps in the Housing



Instrument Implant kit



Manual Screwdriver





<https://alodenti.com/>

09927311262

<https://instagram.com/alodenti.shop>

03134341803

**Instrument**

**Product**

**Photo**

Long Screwdriver  
Ratchet



Short Screwdriver  
Ratchet



Long Screwdriver  
Ratchet



Short Screwdriver  
Ratchet



Screwdriver  
Long Lab



Screwdriver for angle



Screwdriver Locator Abutment  
driver for square



Ratchet Wrench



Universal Ratchet with 4 mm  
connector ( Bone Expander )

